



Color Options for Item: 5REF2030

Note: The PMS colors listed are only intended to provide a general approximation of the tissue color and are not exact matches.

Green on
Green
(GNG)
PMS 7723



Red on
Red
(RNR)
PMS 1795



Silver on
Blue
(SNB)
PMS 2935/
PMS 877



Red
Reflections
(RED)
White/
PMS 032



Gold
Reflections
(GLD)
White/
PMS 871



Silver
Reflections
(SLV)
White/
PMS 877



Rainbow
Reflections
(RRE)



Gold Hot
Spots
(GLS)



Red Hot
Spots
(RDS)




Silver Hot
Spots
(SLS)



Hot Pink
Hot Spots
(HPS)



Black on Black
Hot Spots
(BBS)



Rainbow Hot
Spots
(RHS)



PLEASE NOTE: Colors shown are simulated representations only. Due to the differences in monitors and preferences, there may be slight differences in hue, saturation and tint.